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

Digital clinical thermometer

GB

Thank you for choosing this instrument from TFA.

1. Before you use it

- Please be sure to read the instruction manual carefully.
- Following the instruction manual for use will prevent damage to the device and loss of your statutory rights arising from defects due to incorrect use.
- We shall not be liable for any damage occurring as a result of not following these instructions. Likewise, we take no responsibility for any incorrect readings and for any consequences which may result from them.
- Please take particular note of the safety advice!
- Please look after this manual for future reference.

2. For your safety

- The unit is designed for home use only and suitable exclusively for measuring human body temperature. The product should only be used as described within this instruction.
- Before each measurement check intactness of tip and display. Use the unit only when intact.
- Unauthorised repairs, modifications or changes to the product are prohibited.
- Operation in strong electromagnetic fields, such as next to a mobile phone, can lead to malfunction.
- The minimum measurement time until the beep must be maintained without exception.
- Discuss the measured values with your doctor.

**Caution!**
Risk of injury:

- The instrument contains small parts that could be swallowed by children. Keep this instrument and the battery out of reach of children. Use under parental supervision.
- Batteries and small pieces can be fatal if swallowed. If a battery or a small piece has been swallowed, get medical assistance immediately.
- Batteries must not be thrown into the fire, short-circuited, taken apart or recharged. Risk of explosion!
- Batteries contain harmful acids. Low batteries should be changed as soon as possible to prevent damage caused by a leaking battery.
- Wear chemical-resistant protective gloves and glasses when handling leaked batteries.

**Important information on product safety!**

- Do not place the unit near extreme temperatures, vibration or shock. Do not bend the instrument.
- The thermometer is not watertight. Do not immerse the instrument in water.

3. Elements and buttons**A:** Battery compartment**B:** ON/OFF button**C:** Display for the temperature**D:** Temperature probe**4. Use**

- Disinfect the tip before each use!
- Press the ON/OFF button for switching on and off.
- When switching on a brief beep tone will sound and all LCD segments will light up for a short moment. The instrument starts a self-test.
- The display shows for about 2 seconds the last measured temperature.
- L stays in the display and °C flashes. The unit is ready for use.
- The actually measured temperature is permanently shown on the display and the symbol °C flashes. The measurement is completed when the temperature rises less than 0.1 °C within 16 seconds.
- If the maximum temperature remains stable, an alarm will sound four times and °C is no longer flashing. By extending the duration of measurement after the acoustic signal better measurement results are achieved. For axillary measurements, the recommended measuring time has to be observed despite the beep tone.
- Remove the thermometer carefully. Press the ON/OFF button to switch off the instrument. The instrument will turn off automatically after approx. 10 minutes.

4.1 Measurement methods**Rectal measurement:**

- This is the safest method of measurement. It is especially good for babies and infants. Gently insert the thermometer tip 1-2 cm into the anus. Minimum measuring time: 1 minute.

Oral measurement:

- Place the temperature tip in one of the two pockets under the tongue, left or right at the root of the tongue. Minimum measuring time: 1 minute, mean deviation to the rectal method: 0.3 °C to 0.8 °C.

Axillary measurement (under the arm):

- From the medical point of view this method always provides inaccurate values. Minimum measurement time: 10 minutes, mean deviation to the rectal method: 0,5 °C to 1,5 °C

5. Care and maintenance

- Clean the instrument after each use with a soft damp cloth. For disinfection please use one of the below recommended substances. Any other cleansing substances or cleansing method can cause dysfunction or damage to the instrument. Product names: Edisonite, Gigasept FF, Ethyl alcohol 95,0 %.
- Remove the battery if you do not use the product for a lengthy period.
- Keep the instrument in its case in a dry place.

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GB

5.1 Battery replacement

- Change the battery when the battery symbol  appears on the right hand corner on the display.
- Remove the blue cover and slide the battery compartment forward.
- Remove the battery carefully with a sharp object and insert a new battery LR 41. Make sure the polarity is correct (+ pole above).
- Slide back the battery compartment and replace the blue cover.

6. Troubleshooting

Problem	Solution
No display	→ Switch on the instrument → Ensure battery polarity is correct (+ pole above) → Change battery
Incorrect display	→ Change battery → Check the correct position of the thermometer (see: Measurement methods)
Display L °C	→ Surrounding temperature lower than 32 °C
Display H °C	→ Surrounding temperature higher than 42 °C

If your device fails to work despite these measures contact the supplier from whom you purchased it.

7. Waste disposal

This product has been manufactured using high-grade materials and components which can be recycled and reused.

Never dispose empty batteries and rechargeable batteries in household waste.

As a consumer, you are legally required to take them to your retail store or to appropriate collection sites depending to national or local regulations in order to protect the environment.

The symbols for the heavy metals contained are: Cd=cadmium, Pb=lead, Hg=mercury

This instrument is labelled in accordance with the EU Waste Electrical and Electronic Equipment Directive (WEEE). Please do not dispose of this instrument in household waste. The user is obligated to take end-of-life devices to a designated collection point for the disposal of electrical and electronic equipment, in order to ensure environmentally-compatible disposal.

8. Specifications

Type:	Maximum thermometer
Measuring range:	32 °C...42 °C
L °C	Surrounding temperature lower than 32 °C
H °C	Surrounding temperature higher than 42 °C
Accuracy:	±0,1 °C to 35.5 °C...42.0 °C ±0,2 °C to 32.0 °C...35.5 °C
Resolution:	0.1 °C
Storage range:	-25 °C...+55 °C / humidity: ≤ 95 % rH
Ambient temperature during use:	5 °C...35 °C / humidity: ≤ 80 % rH
Battery:	LR 41 included
Housing dimension:	127 x 10 mm
Weight:	10 g (incl. battery)

Device complies with protection class B.

This instrument complies with the directive for medical devices (93/42/EWG).

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EU Representative: Shanghai International Holding Corp. GmbH (Europe), Eiffestrasse 80, 20537 Hamburg, Germany

Manufacturer: Hangzhou Hua'an Medical & Health Instruments Co., Ltd, Building 2, Baimiao Ind. Park, Wuchang Economic Develop Zone, Hangzhou, 310023, China

No part of this manual may be reproduced without written consent of TFA Dostmann. The technical data are correct at the time of going to print and may change without prior notice.

www.tfa-dostmann.de

05/12

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TFA Dostmann GmbH & Co. KG, Zum Ottersberg 12, D-97877 Wertheim.